

ACI-LeadFree Soldering Summit


SnPb Solder Analysis by X-ray fluorescence
Spectroscopy

Michael Haller
Fischer Technology, Inc.

July 12-13th Philadelphia/PA


The logo for Fischer Technology, Inc. features the word "Fischer" in a bold, white, sans-serif font. A horizontal line is positioned above the letters "i", "s", and "c". A registered trademark symbol (®) is located to the right of the word. The logo is set against a black background with orange bars on either side.

X-Ray fluorescence process



- Incident X-ray Beam strikes sample
- Excitation of characteristic X-ray Fluorescence (photoelectric effect)
- Element Specific characteristic radiation is detected
- Software evaluates Spectrum
- Non-destructive and quick

Hardware Configuration




- 1. X-ray tube
- 2. Cathode
- 3. Anode
- 4. Primary Beam
- 5. Camera Mirror
- 6. Camera Optic
- 7. Video Camer
- 8. Collimator
- 9. Sample Top Layer
- 10. Sample Substrate
- 11. Fluorescence Radiation
- 12. Detector
- 13. Countrate
- 14. Spectrum
- 15. Measurement Results

Fischerscope® X-Ray XDAL / XAN



Pb – detection



- 4 μm Sn coating contaminated with 0.4% Pb on FeNi
- Detection limits of 0.01% possible

Market Requirements-Regulatory Driven

■ Lead Free


- **Directive 2002/95/EC** (Final 27 January 2003)
 - **R**estriction of the use of **H**azardous **S**ubstances in electrical and electronic equipment
- **Directive 2002/96/EC** (Final 27 January 2003)
 - **W**aste **E**lectrical and **E**lectronic **E**quipment
 - Prohibits certain materials in landfills, forces recycling
- **Directive 2000/53/EC** (Final 18 September 2000)
 - **E**nd of **L**ife **V**ehicles
 - Mostly affects Automobiles, but includes language on automotive electronics

Market Requirements - Reliability Driven


- Avoiding Lead Free - Minimum of 3 % Pb to reduce risk of Sn – Whiskering
- End users demand “not pure tin finishes”
 - Specifications
 - Boeing BQA-96-03
 - Requirements for Soldered/Plated Electrical, Electronic Assemblies/Harness/Cables/Components and Mechanical Items
 - Defense Supply Center Columbus
 - MIL-PRF-38534
 - » General performance requirements for hybrid microcircuits, Multi-Chip Modules (MCM)

Typical Sample



- Multi-Layered Chip Capacitor
 - Various Solder alloys
 - Ni or Cu – Barrier Layers
 - Various metallization alloys
 - Ceramic material with many doping elements



- SnPb Solder Terminations
- Ni-Barrier Layer on Terminations
- Ag-Pb-metallization layers
- Ceramic materials



Typical Sample - MLCC's




- MLCC 'S (MultiLayerCeramicCapacitors)
 - Ceramic Base Material
 - Contains many elements Pb,BiAg.
 - Ag in Solder and Paste/Frit
- Application needs to be treated as a layered system
- Any other setup will yield wrong results with potential wrong classification of parts
 - Pure Sn- treated as containing > 3% Pb, Bi, Ag

High Reliability Industry View

- Mil/Aero not required to change by RoHS, but very much affected due to component termination change to 100% Tin
- Tin Whiskers are a real threat to high reliability electronics
- Mitigation of 100% Tin after it sneaks in is far too expensive. Costs will increase, need to be passed on to customers
- Can keep SnPb by working with suppliers and by adding the right equipment
- Flow down requirements to manage 100% Tin are coming soon.
- Can be turned into a competitive advantage

XRF Spectrum of a MLCC



Measurement Results

Application	SnPb[μ"]	Sn [%]	Pb [%]	Ni [μ"]	AgPb[μ"]	Ag [%]	Pb [%]
SnPb/Ni/Ag/Ceramic	97.1 ± 1.52	94.5 ± 0.18	5.48 ± 0.18	68.2 ± 1.1	1313 ± 11.4		
SnPb/Ni/AgPb/Ceramic	74.0 ± 1.43	100.00 ± 0.02	-0.02 ± 0.02	53.5 ± 0.34	1300 ± 18.5	98.3 ± 0.05	1.74 ± 0.05
		Sn I [%]	Pb I [%]	Ni I [%]		Ag I [%]	
SnPbNiAg		12.8 ± 0.17	1.61 ± 0.05	22.3 ± 0.07		63.3 ± 0.20	

■ 3 Different Applications

- Ignoring Pb in Ag-paste/Ceramic
- Define Pb in Ag-paste/Ceramic
- Ignore layered structure of MLCC- measured as bulk

Measurement Results

- Importance of taking layered structure into account
- Significant differences for different applications
- Pure Sn finishes could pass as containing $> 3\%$ Pb, Ag etc
- Bulk analysis not possible for layered samples

Automatic Product Search

Search automatically for the fitting product


Product class: MLCC:Ad

Name	mq (b...	No.	Directory	Date	Measuring mode
SnPb/Ni/AgPb/Ceramic	1.518	5	Fischer	6/23/2...	dcddcc
SnPb/Ni/AgPdPb/Ceramic	1.528	2	Fischer	6/23/2...	dcddccc
SnPb/Ni/AgPd/Ceramic	1.632	1	Fischer	6/20/2...	dcddcc
SnPb/Ni/Ag/Ceramic	1.643	3	Fischer	6/22/2...	dcdd
SnPbNiAg	4.127	4*	Fischer	6/22/2...	Cccc
Au/Pd/Ni/CuFe	20.115	6	Fischer	6/22/2...	dddCc

Measure and search 60 Meas. time Accept Cancel Help

Elemental Mapping of PCB-Board

- Scanning Capabilities with programmable XY-Stage
- Determination of Pb-presence
 - Relevant for RoHS-compliance



Conclusion

- Huge concern in high reliability industry - looking for a solution
- XRF provides a solution
 - Non-destructive
 - No sample preparation
 - Quick (30-60 s)
 - Scanning capabilities
 - Automated Product identification
- XRF requires layered sample measurement capabilities