

## Environmental Chambers

---

### Features

- Stainless steel internal liners and external shells
- Temperature Range:
  - - $-20^{\circ}\text{C}$  to  $177^{\circ}\text{C}$  at sensor after stabilization
- Humidity Range:
  - - 20 to 95% as shown on link below,  $\pm 5\%$  RH at sensor after stabilization
- Includes standard type K thermocouple
- Double plenum feature
- Standard power: 120 VAC, 1ph., 50/60Hz. or 240 VAC, 1ph., 50/60Hz.
- Wide selection of chamber sizes, custom models available



### Accessories

- View port
- Internal light
- Slide assembly
- Load train ports
- Door switch
- Additional options available, contact ATS for more information!