

HSW High Efficiency LED Illumination Handle

- Replacement for conventional light sources
- Lightweight and compact
- LI-Ion Battery / low self-discharge
- Network independent operation
- High grade housing aluminium anodized
- Supervised battery condition
- Adjustable light intensity
- Low heat generation of heat
- Very high light performance
- Simple and rugged mounting to R and PR HSW Endoscopes
- Also lockable to RL Endoscopes and scopes from other manufacturers if RIWO or ACM connection system is used.
- Automatic charger unit enclosed in delivery
- No twisting of the handle in relation to the scope



Technical datas:

Light output:	adjustable in 5 steps from 41000 till 185000 Lux
Operating time:	8h (step 1) 1,2h (step 5)
Max. power :	5W
Colour temperature:	5500°K
Durability LED:	50 000h
Charge time:	2,5 h