

COMMERCIAL - BENCHTOP ULTRASONIC CLEANERS



The FXP series of ultrasonic cleaners are designed for Medical, Dental, Laboratory and where all forms of contamination are lightly deposited.

All units incorporate Lateral Oscillation Technology for uniform cleaning results. The internal basket system allows the lid to seal the working solution avoiding evaporation and reducing ultrasonic acoustic levels. The clear and easy design with compact and quiet operation make these products the market leaders.

PEARLS, AMBER AND OPALS ARE POROUS GEMS AND NOT SUITABLE FOR THE ULTRASONIC CLEANER. THIS TYPE OF JEWELLERY SHOULD ONLY BE CLEANED USING A SOFT CLOTH WITH A MILD SOLUTION.



FXP4 **JEWELLERY CLEANER**

Chamber capacity: 1.1L.
Operating Liquid capacity 0.75 Litres
Glass Cleaning Platform;
Length 130mm
Width 130mm
Height 40mm



FXP4M

Chamber Volume 1.1 Litres
Operating Liquid capacity 0.75 Litres
Basket Size;
Length 120mm
Width 100mm
Height 40mm



FXP8M

Chamber Volume 1.7 Litres
Operating Liquid capacity 1.2 Litres
Basket Size;
Length 120mm
Width 100mm
Height 65mm

Ultrasonic Cleaning Tanks

FSS High Performance - ST Industrial - SIEVE Cleaners



FXP8D

Chamber Volume 1.7 Litres
Operating Liquid capacity 1.2 Litres
Basket Size;
Length 120mm
Width 100mm
Height 65mm



FXP10M

Chamber Volume 2.7 Litres
Operating Liquid capacity 2 Litres
Basket Size;
Length 210mm
Width 100mm
Height 65mm



FXP10D

Chamber Volume 2.7 Litres
Operating Liquid capacity 2 Litres
Basket Size;
Length 210mm
Width 100mm
Height 65mm



FXP10MH

Chamber Volume 2.7 Litres
Operating Liquid capacity 2 Litres
Basket Size;
Length 210mm
Width 100mm
Height 65mm



FXP10DH

Chamber Volume 2.7 Litres
Operating Liquid capacity 2 Litres
Basket Size;
Length 210mm
Width 100mm
Height 65mm



FXP12M

Chamber Volume 5.6 Litres
Operating Liquid capacity 4.5 Litres

Basket Size;

Length 250mm
Width 110mm
Height 100mm



FXP12D

Chamber Volume 5.6 Litres
Operating Liquid capacity 4.5 Litres

Basket Size;

Length 250mm
Width 110mm
Height 100mm



FXP12MH

Chamber Volume 5.6 Litres
Operating Liquid capacity 4.5 Litres

Basket Size;

Length 250mm
Width 110mm
Height 100mm



FXP12DH

Chamber Volume 5.6 Litres
Operating Liquid capacity 4.5 Litres

Basket Size;

Length 250mm
Width 110mm
Height 100mm



FXP14M

Chamber Volume 10 Litres
Operating Liquid capacity 8 Litres

Basket Size;

Length 250mm
Width 200mm
Height 100mm



FXP14D

Chamber Volume 10 Litres
Operating Liquid capacity 8 Litres
Basket Size;
Length 250mm
Width 200mm
Height 100mm



FXP14MH

Chamber Volume 10 Litres
Operating Liquid capacity 8 Litres
Basket Size;
Length 250mm
Width 200mm
Height 100mm



FXP14DH

Chamber Volume 10 Litres
Operating Liquid capacity 8 Litres
Basket Size;
Length 250mm
Width 200mm
Height 100mm



FXP16M

Chamber Volume 14 Litres
Operating Liquid capacity 11 Litres
Basket Size;
Length 280mm
Width 260mm
Height 100mm



FXP16D

Chamber Volume 14 Litres
Operating Liquid capacity 11 Litres
Basket Size;
Length 280mm
Width 260mm
Height 110mm



FXP16MH

Chamber Volume 14 Litres
Operating Liquid capacity 11 Litres
Basket Size;
Length 280mm
Width 260mm
Height 100mm



FXP16DH

Chamber Volume 14 Litres
Operating Liquid capacity 11 Litres
Basket Size;
Length 280mm
Width 260mm
Height 100mm



FXP20M

Chamber Volume 21.5 Litres
Operating Liquid capacity 18 Litres
Basket Size;
Length 450mm
Width 250mm
Height 100mm



FXP20D

Chamber Volume 21.5 Litres
Operating Liquid capacity 18 Litres
Basket Size;
Length 450mm
Width 250mm
Height 100mm



FXP20MH

Chamber Volume 21.5 Litres
Operating Liquid capacity 18 Litres
Basket Size;
Length 450mm
Width 250mm
Height 100mm



FXP20DH

Chamber Volume 21.5 Litres
Operating Liquid capacity 18 Litres
Basket Size;
Length 450mm
Width 250mm
Height 100mm




FXP40M

Chamber Volume 61 Litres
Operating Liquid capacity 52 Litres
Basket Size;
Length 550mm
Width 450mm
Height 140mm




FXP40D

Chamber Volume 61 Litres
Operating Liquid capacity 52 Litres
Basket Size;
Length 550mm
Width 450mm
Height 140mm



FXP40MH

Chamber Volume 61 Litres
Operating Liquid capacity 52 Litres
Basket Size;
Length 550mm
Width 450mm
Height 140mm



FXP40DH

Chamber Volume 61 Litres
Operating Liquid capacity 52 Litres
Basket Size;
Length 550mm
Width 450mm
Height 140mm